




CA194 rabbit pAb

Catalog No	YP-Ab-08481
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	C1orf194
Immunogen	Synthesized peptide derived from human CA194 AA range: 55-105
Specificity	This antibody detects endogenous levels of CA194 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CA194;C1orf194
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Cytoplasm, cytoskeleton .
Tissue Specificity	Expressed in cerebrum, cerebellum, gastrocnemius muscle, spinal cord and lung tissues.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from DU145 cells,
primary antibody was diluted at 1:1000, 4° over night