




CE025 rabbit pAb

Catalog No	YP-Ab-08488
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	C5orf25
Immunogen	Synthesized peptide derived from human CE025 AA range: 369-419
Specificity	This antibody detects endogenous levels of CE025 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CE025;C5orf25
Observed Band	
Calculated Molecular Weight	
Cell Pathway	[Isoform 5]: Cytoplasm . Cytoplasm, myofibril, sarcomere .
Tissue Specificity	[Isoform 5]: Skeletal muscle.
Function	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from 3T3 cells, primary antibody was diluted at 1:1000, 4° over night