



TBG2 rabbit pAb


Catalog No	YP-Ab-08498
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	TUBG2
Immunogen	Synthesized peptide derived from human TBG2 AA range: 113-163
Specificity	This antibody detects endogenous levels of TBG2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Gamma 1 tubulin;Gamma 2 tubulin;Gamma Tubulin 1;Gamma Tubulin 2; Gamma tubulin complex component 1;Gamma-2-tubulin;GCP 1;GCP1; MGC131994;TBG2_HUMAN;TUBG;TUBG1;TUBG2;TUBGCP1;Tubulin gamma 1 chain;Tubulin gamma 2 chain;Tubulin gamma complex-associated protein 1; Tubulin gamma-2 chain;tubulin;gamma 1;gamma 2;gamma polypeptide;Xgam; TBG2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm, cytoskeleton, microtubule organizing center, centrosome .
Tissue Specificity	
Function	function:Tubulin is the major constituent of microtubules. Gamma tubulin is found at microtubule organizing centers (MTOC) such as the spindle poles or the centrosome, suggesting that it is involved in the minus-end nucleation of microtubule assembly.,similarity:Belongs to the tubulin family.,
Background	

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4°over night