



UBR3 rabbit pAb

Catalog No	YP-Ab-08500
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	UBR3 KIAA2024 ZNF650
Immunogen	Synthesized peptide derived from human UBR3 AA range: 1591-1641
Specificity	This antibody detects endogenous levels of UBR3 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	bA49A4.1;C6orf133;dJ242G1.1;dJ392M17.3;DKFZp686C8114;E3 ubiquitin protein ligase UBR2;KIAA349;MGC71112;N recognin 2;OTTHUMP1643; OTTHUMP39767;OTTHUMP39768;RP3 392M17.3;Ubiquitin ligase E3 alpha II; Ubiquitin protein ligase E3 alpha 2;Ubiquitin protein ligase E3 alpha II;Ubiquitin protein ligase E3 component n recognin 2;UBR 2;UBR3;UBR3 KIAA2024 ZNF650
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	function:E3 ubiquitin-protein ligase which is a component of the N-end rule pathway. Does not bind to proteins bearing specific N-terminal residues that are destabilizing according to the N-end rule, leading to their ubiquitination and subsequent degradation.,pathway:Protein modification; protein ubiquitination., sequence caution:Translated as Gln.,similarity:Belongs to the UBR1 family., similarity:Contains 1 RING-type zinc finger.,similarity:Contains 1 UBR-type zinc finger.,subunit:Interacts with UBE2A and UBE2B.,

Background


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from THP-1 cells, primary antibody was diluted at 1:1000, 4° over night