




COLL1 rabbit pAb

Catalog No	YP-Ab-08502
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	CLPSL1 C6orf127
Immunogen	Synthesized peptide derived from human COLL1 AA range: 35-85
Specificity	This antibody detects endogenous levels of COLL1 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	coll 1;Hsema I;Hsema III;Sema III;SEMA1;SEMA3A;SEMAD;SEMAIII;SEMAL; Semaphorin 3A;Semaphorin III;SemD;COLL1;CLPSL1 C6orf127
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted .
Tissue Specificity	Exclusively expressed in epididymis, in the corpus region.
Function	similarity:Belongs to the colipase family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from A431 cells,
primary antibody was diluted at 1:1000, 4° over night