




CDR2L rabbit pAb

Catalog No	YP-Ab-08508
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	CDR2L HUMPPA
Immunogen	Synthesized peptide derived from human CDR2L AA range: 319-369
Specificity	This antibody detects endogenous levels of CDR2L at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Paraneoplastic 62 kDa antigen;HUMPPA;Cerebellar degeneration-related protein 2-like;CDR2L;CDR2L HUMPPA
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	similarity:Belongs to the CDR2 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from U2OS cells,
primary antibody was diluted at 1:1000, 4° over night