




MYL10 rabbit pAb

Catalog No	YP-Ab-08514
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	MYL10 MYLC2PL PLRLC
Immunogen	Synthesized peptide derived from human MYL10 AA range: 78-128
Specificity	This antibody detects endogenous levels of MYL10 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MYLC2PL;Myosin light chain 2;lymphocyte-specific;PLRLC;Precursor lymphocyte-specific regulatory light chain;MYL10;MYL10 MYLC2PL PLRLC
Observed Band	
Calculated Molecular Weight	
Cell Pathway	mitochondrion,cytosol,
Tissue Specificity	
Function	miscellaneous:This chain binds calcium.,similarity:Contains 3 EF-hand domains., subunit:Myosin is an hexamer of 2 heavy chains and 4 light chains.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HpeG2 cells,
primary antibody was diluted at 1:1000, 4° over night