



S6OS1 rabbit pAb


Catalog No	YP-Ab-08522
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	SIX6OS1 C14orf39
Immunogen	Synthesized peptide derived from human S6OS1 AA range: 511-561
Specificity	This antibody detects endogenous levels of S6OS1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C14orf39;Chromosome 14 open reading frame 39;Protein SIX6OS1;S6OS1_HUMAN;Six6 opposite strand transcript 1;Six6os1;S6OS1;SIX6OS1 C14orf39
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Chromosome . Component of the central element of the synaptonemal complex. In spermatocytes, detected from zygonema to pachynema and localizes along synapsed lateral elements. Loading to the central element of the synaptonemal complex is dependent on the assembly of the tripartite synaptonemal complex structure that occurs upon synapsis between homologous chromosomes. .
Tissue Specificity	Highest expression in retina, skeletal muscle, testis and colon.
Function	function:May be involved in eye development through regulation of the SIX6 transcription factor.,tissue specificity:Highest expression in retina, skeletal muscle and colon.,
Background	

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from SH-SY5Y cells, primary antibody was diluted at 1:1000, 4° over night