




CALCB rabbit pAb

Catalog No	YP-Ab-08526
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	CALCB CALC2
Immunogen	Synthesized peptide derived from human CALCB AA range: 45-95
Specificity	This antibody detects endogenous levels of CALCB at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Beta type CGRP;CALC2;CGRP II;CGRP2;CALCB;CALCB CALC2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted.
Tissue Specificity	
Function	function:CGRP induces vasodilation. It dilates a variety of vessels including the coronary, cerebral and systemic vasculature. Its abundance in the CNS also points toward a neurotransmitter or neuromodulator role.,similarity:Belongs to the calcitonin family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from A549 cells,
primary antibody was diluted at 1:1000, 4°over night