




SPRR3 rabbit pAb

Catalog No	YP-Ab-08528
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	SPRR3 SPRC
Immunogen	Synthesized peptide derived from human SPRR3 AA range: 23-73
Specificity	This antibody detects endogenous levels of SPRR3 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SPRC;Small proline-rich protein 3;SPRR3;SPRR3 SPRC
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm.
Tissue Specificity	
Function	function:Cross-linked envelope protein of keratinocytes.,induction:During squamous differentiation of epidermal keratinocytes.,similarity:Belongs to the cornifin (SPRR) family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from Hela cells,
primary antibody was diluted at 1:1000, 4° over night