



ZN350 rabbit pAb


Catalog No	YP-Ab-08531
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	ZNF350 ZBRK1
Immunogen	Synthesized peptide derived from human ZN350 AA range: 346-396
Specificity	This antibody detects endogenous levels of ZN350 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KRAB zinc finger protein ZFQR;ZBRK1;ZFQR;Zinc finger and BRCA1 interacting protein with a KRAB domain 1;Zinc finger and BRCA1-interacting protein with a KRAB domain 1;Zinc finger protein 35;Zinc finger protein ZBRK1; ZN35_HUMAN;ZNF35;ZN350;ZNF350 ZBRK1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus. Nucleus matrix. Associated with the nuclear matrix.
Tissue Specificity	Widely expressed.
Function	function:Transcriptional repressor. Binds to a specific sequence, 5'-GGGxxxCAGxxxTTT-3', within GADD45 intron 3.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 KRAB domain.,similarity:Contains 8 C2H2-type zinc fingers.,subcellular location:Associated with the nuclear matrix.,subunit:Interacts with BRCA1.,tissue specificity:Widely expressed.,
Background	

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products ImagesWestern blot analysis of lysates from U2OS cells,
primary antibody was diluted at 1:1000, 4° over night