




LLPH rabbit pAb

Catalog No	YP-Ab-08535
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	LLPH C12orf31
Immunogen	Synthesized peptide derived from human LLPH AA range: 18-68
Specificity	This antibody detects endogenous levels of LLPH at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	hLLP;C12orf31;LLPH;LLPH C12orf31
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus, nucleolus . Chromosome . Cell-permeable protein. 22 hours after injection in the hippocampal area CA1, internalized by most cells at the injection site (By similarity). Localizes at the chromosome periphery during mitosis (PubMed:20813266). .
Tissue Specificity	
Function	similarity:Belongs to the learning-associated protein family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MCF-7 cells,
primary antibody was diluted at 1:1000, 4°over night