



FRG2 rabbit pAb


Catalog No	YP-Ab-08537
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	FRG2 FRG2A
Immunogen	Synthesized peptide derived from human FRG2 AA range: 104-154
Specificity	This antibody detects endogenous levels of FRG2 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FARP2;FARP2_HUMAN;FERM;FERM domain including RhoGEF;RhoGEF and pleckstrin domain protein 2;RhoGEF and pleckstrin domain-containing protein 2; FGD1 related Cdc42 GEF;FIR;FRG;KIAA793;PH domain-containing family C member 3;Pleckstrin homology domain-containing family C member 3; PLEKHC3;FRG2;FRG2 FRG2A
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	Expression is undetectable in all tissues tested except for differentiating myoblasts of FSHD patients, which display low, yet distinct levels of expression, partly from FRG2, but predominantly originating from its homolog on chromosome 10.
Function	similarity:Belongs to the FRG2 family.,tissue specificity:Expression is undetectable in all tissues tested except for differentiating myoblasts of FSHD patients, which display low, yet distinct levels of expression, partly from FRG2, but predominantly originating from its homolog on chromosome 10.,
Background	

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from HCT116 cells, primary antibody was diluted at 1:1000, 4° over night