



TULP4 rabbit pAb


Catalog No	YP-Ab-08542
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	TULP4 KIAA1397 TUBL4 TUSP
Immunogen	Synthesized peptide derived from human TULP4 AA range: 798-848
Specificity	This antibody detects endogenous levels of TULP4 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TUSP;TUBL4;Tubby-related protein 4;Tubby-like protein 4;KIAA1397;TULP4;TULP4 KIAA1397 TUBL4 TUSP
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm.
Tissue Specificity	Expressed mainly in the brain, skeletal muscle, testis and kidney.
Function	domain:The SOCS box domain mediates the interaction with the Elongin BC complex, an adapter module in different E3 ubiquitin ligase complexes., function:May be a substrate-recognition component of a SCF-like ECS (Elongin-Cullin-SOCS-box protein) E3 ubiquitin ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins., pathway:Protein modification; protein ubiquitination.,similarity:Belongs to the TUB family.,similarity:Contains 1 SOCS box domain.,similarity:Contains 3 WD repeats.,tissue specificity:Expressed mainly in the brain, skeletal muscle, testis and kidney.,
Background	

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from MDA-MB cells, primary antibody was diluted at 1:1000, 4° over night