




CC170 rabbit pAb

Catalog No	YP-Ab-08543
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	CCDC170 C6orf97
Immunogen	Synthesized peptide derived from human CC170 AA range: 68-118
Specificity	This antibody detects endogenous levels of CC170 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	bA282P11.1;CC17_HUMAN;CCDC17;Chromosome 6 open reading frame 97; Coiled coil domain containing protein C6orf97;Coiled-coil domain-containing protein 17;FLJ2335;Hypothetical protein LOC8129;LOC8129;CC170;CCDC170 C6orf97
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Golgi apparatus .
Tissue Specificity	
Function	
Background	The function of this gene and its encoded protein is not known. Several genome-wide association studies have implicated the region around this gene to be involved in breast cancer and bone mineral density, but no link to this specific gene has been found. [provided by RefSeq, May 2010],
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from K562 cells,
primary antibody was diluted at 1:1000, 4°over night