




RISC rabbit pAb

Catalog No	YP-Ab-08550
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	SCPEP1 RISC SCP1 MSTP034 UNQ265/PRO302
Immunogen	Synthesized peptide derived from human RISC AA range: 151-201
Specificity	This antibody detects endogenous levels of RISC at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HSCP1; RISC; SCP1; SCPEP1; serine carboxypeptidase 1; UNQ265/PRO302; SCPEP1 RISC SCP1 MSTP034 UNQ265/PRO302
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted .
Tissue Specificity	
Function	function: May be involved in vascular wall and kidney homeostasis., similarity: Belongs to the peptidase S10 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night