




ACOT6 rabbit pAb

Catalog No	YP-Ab-08551
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	ACOT6 C14orf42
Immunogen	Synthesized peptide derived from human ACOT6 AA range: 116-166
Specificity	This antibody detects endogenous levels of ACOT6 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ACOT6;acyl CoA thioesterase 6;c14_5530;C14orf42;ACOT6 C14orf42
Observed Band	
Calculated Molecular Weight	
Cell Pathway	[Isoform 1]: Peroxisome . Localization to the peroxisome is uncertain since the potential C-terminal peroxisome targeting signal found in the mouse ortholog is not perfectly conserved. .; [Isoform 2]: Cytoplasm . Recombinant N-terminally GFP-tagged protein localizes to the cytosol. .
Tissue Specificity	
Function	caution:May lack catalytic activity, even though the active site residues are conserved.,function:Unknown.,similarity:Belongs to the C/M/P thioester hydrolase family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from AD293 cells, primary antibody was diluted at 1:1000, 4° over night