



SSX6 rabbit pAb

Catalog No	YP-Ab-08552
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	SSX6
Immunogen	Synthesized peptide derived from human SSX6 AA range: 119-169
Specificity	This antibody detects endogenous levels of SSX6 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	dJ54B2.1;OTTHUMP23236;OTTHUMP61616;Protein SSX6;psiSSX2;RP11-38O23.6;SSX 6;SSXP2;Synovial sarcoma X breakpoint 6;Synovial sarcoma X breakpoint 6 protein;SSX6
Observed Band	
Calculated Molecular Weight	
Cell Pathway	nucleus,
Tissue Specificity	Not detected in any normal tissues. Expressed in a melanoma cell line.
Function	function:Could act as a modulator of transcription.,similarity:Belongs to the SSX family.,similarity:Contains 1 KRAB-related domain.,tissue specificity:Not detected in any normal tissues. Expressed in a melanoma cell line.,
Background	This gene belongs to the family of highly homologous synovial sarcoma X (SSX) breakpoint proteins. These proteins may function as transcriptional repressors. They are also capable of eliciting spontaneously humoral and cellular immune responses in cancer patients, and are potentially useful targets in cancer vaccine-based immunotherapy. SSX1, SSX2 and SSX4 genes have been involved in the t(X;18) translocation characteristically found in all synovial sarcomas. This gene is classified as a pseudogene because a splice donor in



the 3' UTR has changed compared to other family members, rendering the transcript a candidate for nonsense-mediated mRNA decay (NMD). [provided by RefSeq, Aug 2009],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from HeLa cells, primary antibody was diluted at 1:1000, 4° over night