




DUS1L rabbit pAb

Catalog No	YP-Ab-08555
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	DUS1L
Immunogen	Synthesized peptide derived from human DUS1L AA range: 262-312
Specificity	This antibody detects endogenous levels of DUS1L at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DUS1; DUS1L; PP3111
Observed Band	
Calculated Molecular Weight	
Cell Pathway	cytosol,
Tissue Specificity	
Function	cofactor: FAD., function: Catalyzes the synthesis of dihydrouridine, a modified base found in the D-loop of most tRNAs., similarity: Belongs to the dus family. Dus1 subfamily.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from A431 cells,
primary antibody was diluted at 1:1000, 4° over night