



ACHA2 rabbit pAb

Catalog No	YP-Ab-08558
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	CHRNA2
Immunogen	Synthesized peptide derived from human ACHA2 AA range: 219-269
Specificity	This antibody detects endogenous levels of ACHA2 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ACHA2; CHRNA2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell junction, synapse, postsynaptic cell membrane; Multi-pass membrane protein. Cell membrane; Multi-pass membrane protein.
Tissue Specificity	
Function	disease: Defects in CHRNA2 are the cause of nocturnal frontal lobe epilepsy type 4 (ENFL4) [MIM:610353]. ENFL4 is an autosomal dominant epilepsy characterized by nocturnal seizures associated with fear sensation, tongue movements, and nocturnal wandering, closely resembling nightmares and sleep walking. function: After binding acetylcholine, the AChR responds by an extensive change in conformation that affects all subunits and leads to opening of an ion-conducting channel across the plasma membrane. similarity: Belongs to the ligand-gated ionic channel (TC 1.A.9) family. subunit: Neuronal AChR seems to be composed of two different types of subunits: alpha and non-alpha (beta). Alpha-2 subunit can be combined to beta-2 or beta-4 to give rise to functional receptors.



Background

Nicotinic acetylcholine receptors (nAChRs) are ligand-gated ion channels formed by a pentameric arrangement of alpha and beta subunits to create distinct muscle and neuronal receptors. Neuronal receptors are found throughout the peripheral and central nervous system where they are involved in fast synaptic transmission. This gene encodes an alpha subunit that is widely expressed in the brain. The proposed structure for nAChR subunits is a conserved N-terminal extracellular domain followed by three conserved transmembrane domains, a variable cytoplasmic loop, a fourth conserved transmembrane domain, and a short C-terminal extracellular region. Mutations in this gene cause autosomal dominant nocturnal frontal lobe epilepsy type 4. Single nucleotide polymorphisms (SNPs) in this gene have been associated with nicotine dependence. [provided by RefSeq, Nov 2009],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4° over night