




FA53B rabbit pAb

Catalog No	YP-Ab-08561
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	FAM53B KIAA0140
Immunogen	Synthesized peptide derived from human FA53B AA range: 260-310
Specificity	This antibody detects endogenous levels of FA53B at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FA53B;FAM53B KIAA0140
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	Detected in skeletal muscle, kidney, spleen, thyroid, testis, ovary, small intestine, colon and peripheral blood.
Function	similarity:Belongs to the FAM53 family.,tissue specificity:Detected in skeletal muscle, kidney, spleen, thyroid, testis, ovary, small intestine, colon and peripheral blood.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from THP-1 cells,
primary antibody was diluted at 1:1000, 4° over night