



ENOPH rabbit pAb

Catalog No	YP-Ab-08563
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	ENOPH1 MASA MSTP145
Immunogen	Synthesized peptide derived from human ENOPH AA range: 112-162
Specificity	This antibody detects endogenous levels of ENOPH at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	2 3 diketo 5 methylthio 1 phosphopentane phosphatase;2,3-diketo-5-methylthio-1-phosphopentane phosphatase;Acireductone synthase;DKFZp586M524;E 1 enzyme;E1;Enolase phosphatase 1;Enolase phosphatase E1;Enolase-phosphatase E1;ENOPH HUMAN;Enoph1;FLJ12594;MASA homolog;MST145; MSTP145 protein;ENOPH;ENOPH1 MASA MSTP145
Observed Band	29kD
Calculated Molecular Weight	29kD
Cell Pathway	Cytoplasm . Nucleus .
Tissue Specificity	
Function	catalytic activity:5-(methylthio)-2,3-dioxopentyl phosphate + H(2)O = 1,2-dihydroxy-5-(methylthio)pent-1-en-3-one + phosphate.,cofactor:Binds 1 magnesium ion per subunit.,function:Bifunctional enzyme that enolizes the substrate to form the intermediate 2-hydroxy-5-(methylthio)-3-oxopent-1-enyl phosphate, which is then dephosphorylated to form the acireductone 1,2-dihydroxy-5-(methylthio)pent-1-en-3-one.,pathway:Amino-acid biosynthesis; L-methionine biosynthesis via salvage pathway; L-methionine from (S)-methyl-5-thio-alpha-D-ribose 1-phosphate: step 3/6.,pathway:Amino-acid biosynthesis; L-methionine biosynthesis via salvage pathway; L-methionine from (S)-methyl-5-



thio-alpha-D-ribose 1-phosphate: step 4/6.,similarity:Belongs to the HAD-like hydrolase superfamily. MasA/mtnC family.,subunit:Monomer.,

Background


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from DU145 cells, primary antibody was diluted at 1:1000, 4° over night