




TIKI2 rabbit pAb

Catalog No	YP-Ab-08566
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	TRABD2B TIKI2
Immunogen	Synthesized peptide derived from human TIKI2 AA range: 11-61
Specificity	This antibody detects endogenous levels of TIKI2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LOC388630;Tiki2;TRABD2B;UPF0632 protein A;TRABD2B TIKI2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Single-pass type I membrane protein .
Tissue Specificity	
Function	similarity:Belongs to the UPF0632 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from K562 cells,
primary antibody was diluted at 1:1000, 4° over night