




SLIP rabbit pAb

Catalog No	YP-Ab-08567
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	NUGGC C8orf80 HMFN0672
Immunogen	Synthesized peptide derived from human SLIP AA range: 530-580
Specificity	This antibody detects endogenous levels of SLIP at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SLIP;NUGGC C8orf80 HMFN0672
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus speckle .
Tissue Specificity	Expressed in germinal center B-cell and in lymphomas derived from germinal center B-cell.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from 293T cells,
primary antibody was diluted at 1:1000, 4° over night