




MUC22 rabbit pAb

Catalog No	YP-Ab-08568
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	MUC22 G4 PBMUCL1
Immunogen	Synthesized peptide derived from human MUC22 AA range: 286-336
Specificity	This antibody detects endogenous levels of MUC22 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MUC22;MUC22 G4 PBMUCL1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Single-pass membrane protein .
Tissue Specificity	Expressed in lung by serous cells of the submucosal gland (at protein level). Detected in the placenta, lung and testis.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from U2OS cells,
primary antibody was diluted at 1:1000, 4° over night