




MORN2 rabbit pAb

Catalog No	YP-Ab-08589
Isotype	IgG
Reactivity	Human; Mouse
Applications	IHC;IF
Gene Name	MORN2 MOPT
Immunogen	Synthesized peptide derived from human MORN2 AA range: 25-75
Specificity	This antibody detects endogenous levels of MORN2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:50-200. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BLOCK27;MOPT;MORN motif protein in testis;MORN repeat containing 2; MORN2;MORN2 MOPT
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasmic vesicle, secretory vesicle, acrosome . Nucleus . Associated with the acrosome during post-meiotic stages of spermatogenesis. Appears to migrate to the nucleus during the final step of spermiogenesis. .
Tissue Specificity	
Function	similarity:Contains 2 MORN repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Immunohistochemical analysis of paraffin-embedded human Gastric adenocarcinoma. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).