




RNAS8 rabbit pAb

Catalog No	YP-Ab-08592
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	RNASE8
Immunogen	Synthesized peptide derived from human RNAS8 AA range: 27-77
Specificity	This antibody detects endogenous levels of RNAS8 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RNAS8;RNASE8
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted .
Tissue Specificity	Expressed prominently in the placenta and is not detected in any other tissues examined.
Function	function:Has a low ribonuclease activity.,similarity:Belongs to the pancreatic ribonuclease family.,tissue specificity:Expressed prominently in the placenta and is not detected in any other tissues examined.,
Background	The protein encoded by this gene belongs to the pancreatic ribonuclease family, a subset of the ribonuclease A superfamily. The protein is expressed prominently in the placenta and exhibits antimicrobial activity against pathogenic bacteria and fungi. [provided by RefSeq, Oct 2014],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from MCF-7 cells, primary antibody was diluted at 1:1000, 4° over night