



YN3759

ZG16 rabbit pAb


Catalog No	YP-Ab-08601
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	ZG16
Immunogen	Synthesized peptide derived from human ZG16 AA range: 59-109
Specificity	This antibody detects endogenous levels of ZG16 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FLJ43571;FLJ92276;hZG16;JCLN;JCLN1;ZG16
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted, extracellular space, extracellular matrix . Zymogen granule lumen . Golgi apparatus lumen .
Tissue Specificity	Highly expressed in liver. Detected at lower levels in colon, ileum and jejunum.
Function	function:May act as a linker molecule between the submembranous matrix on the luminal side of zymogen granule membrane (ZGM) and aggregated secretory proteins during granule formation in the TGN.,similarity:Belongs to the jacalin lectin family.,subcellular location:Stored in zymogen granules.,
Background	

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products ImagesWestern blot analysis of lysates from HeLa cells,
primary antibody was diluted at 1:1000, 4°over night