




SFR1 rabbit pAb

Catalog No	YP-Ab-08607
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	SFR1 C10orf78 MEI5 MEIR5
Immunogen	Synthesized peptide derived from human SFR1 AA range: 110-160
Specificity	This antibody detects endogenous levels of SFR1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C10orf78; MEI5; Meiosis protein 5 homolog; MEIR5; SWI5 Dependent Homologous Recombination Repair Protein 1; SFR1; SFR1 C10orf78 MEI5 MEIR5
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus . Colocalizes with ESR1 in the nucleus. .
Tissue Specificity	Widely expressed.
Function	caution: It is uncertain whether Met-1 or Met-76 is the initiator.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from Hela cells,
primary antibody was diluted at 1:1000, 4° over night