



S61A1 rabbit pAb

Catalog No	YP-Ab-08614
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	SEC61A1 SEC61A
Immunogen	Synthesized peptide derived from human S61A1 AA range: 333-383
Specificity	This antibody detects endogenous levels of S61A1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	chunp6898;fb62c11;HSEC61;protein transport protein SEC61 alpha subunit; protein transport protein Sec61 subunit alpha;Protein transport protein Sec61 subunit alpha isoform 1;S61A1_HUMAN;SEC 61;SEC 61A;Sec61 alpha 1; Sec61 alpha 1 subunit;Sec61 alpha 1 subunit (S. cerevisiae);Sec61 alpha-1; SEC61;Sec61 homolog;SEC61A1;sec61aa;S61A1;SEC61A1 SEC61A
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Endoplasmic reticulum membrane ; Multi-pass membrane protein . Localizes exclusively in granular structures in the endoplasmic reticulum (ER). .
Tissue Specificity	Expressed in proximal and distal tubules in kidney (at protein level).
Function	function:Plays a crucial role in the insertion of secretory and membrane polypeptides into the ER. Required for assembly of membrane and secretory proteins. Tightly associated with membrane-bound ribosomes, either directly or through adapter proteins.,similarity:Belongs to the secY/SEC61-alpha family., subunit:Heterotrimeric complex composed of SEC61-alpha, SEC61-beta and SEC61-gamma.,



Background

The protein encoded by this gene belongs to the SECY/SEC61- alpha family. It appears to play a crucial role in the insertion of secretory and membrane polypeptides into the endoplasmic reticulum. This protein found to be tightly associated with membrane-bound ribosomes, either directly or through adaptor proteins. This gene encodes an alpha subunit of the heteromeric SEC61 complex, which also contains beta and gamma subunits. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from DU145 cells, primary antibody was diluted at 1:1000, 4° over night