



CL12B rabbit pAb

Catalog No	YP-Ab-08616
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	CLEC12B UNQ5782/PRO16089
Immunogen	Synthesized peptide derived from human CL12B AA range: 111-161
Specificity	This antibody detects endogenous levels of CL12B at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CL12B;CLEC12B UNQ5782/PRO16089
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Single-pass type II membrane protein .
Tissue Specificity	Detected in colon, heart, kidney, liver, lung, mammary gland, ovary, spleen and testis (PubMed:17562706). Expressed in melanocytes (at protein level) (PubMed:34310951).
Function	domain:Contains 1 copy of a cytoplasmic motif that is referred to as the immunoreceptor tyrosine-based inhibitor motif (ITIM). This motif is involved in modulation of cellular responses. The phosphorylated ITIM motif can bind the SH2 domain of several SH2-containing phosphatases.,function:Cell surface receptor that modulates signaling cascades and mediates tyrosine phosphorylation of target MAP kinases.,function:Cell surface receptor that protects target cells against natural killer cell-mediated lysis. Modulates signaling cascades and mediates tyrosine phosphorylation of target MAP kinases.,induction:Down-regulated in activated leukocytes recruited to a site of inflammation.,PTM:Highly N-glycosylated. Glycosylation varies between cell types.,similarity:Contains 1 C-type lectin domain.,subcellular location:Ligand



binding leads to internalization.,subunit:Homodimer. Interacts with PTPN6 a

Background


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from HeLa cells, primary antibody was diluted at 1:1000, 4°over night