




ANR42 rabbit pAb

Catalog No	YP-Ab-08624
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	ANKRD42
Immunogen	Synthesized peptide derived from human ANR42 AA range: 11-61
Specificity	This antibody detects endogenous levels of ANR42 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ANR42;ANKRD42
Observed Band	
Calculated Molecular Weight	
Cell Pathway	ribosome,
Tissue Specificity	
Function	similarity:Contains 9 ANK repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night