




CT203 rabbit pAb

Catalog No	YP-Ab-08635
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	C20orf203
Immunogen	Synthesized peptide derived from human CT203 AA range: 10-60
Specificity	This antibody detects endogenous levels of CT203 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CT203;C20orf203
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm .
Tissue Specificity	Expressed most abundantly in the brain at protein level. Present in cortex, cerebellum and midbrain. Found in neurons. Elevated expressions detected in Alzheimer brain samples. Also expressed in testis.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from 293T cells,
primary antibody was diluted at 1:1000, 4° over night