



H2AW rabbit pAb

Catalog No	YP-Ab-08642
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	H2AFY2 MACROH2A2
Immunogen	Synthesized peptide derived from human H2AW AA range: 151-201
Specificity	This antibody detects endogenous levels of H2AW at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HIST3H2A;histone;histone cluster 3;H2a;Histone H2A;Histone H2A type 3; H2AW;H2AFY2 MACROH2A2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus . Chromosome . Enriched in inactive X chromosome chromatin (PubMed:11331621, PubMed:11262398) and in senescence-associated heterochromatin (PubMed:15621527). .
Tissue Specificity	
Function	function:Variant histone H2A which replaces conventional H2A in a subset of nucleosomes where it represses transcription. Nucleosomes wrap and compact DNA into chromatin, limiting DNA accessibility to the cellular machineries which require DNA as a template. Histones thereby play a central role in transcription regulation, DNA repair, DNA replication and chromosomal stability. DNA accessibility is regulated via a complex set of post-translational modifications of histones, also called histone code, and nucleosome remodeling. May be involved in stable X chromosome inactivation.,similarity:Contains 1 histone H2A domain.,similarity:Contains 1 Macro domain.,subcellular location:Enriched in inactive X chromosome chromatin and in senescence-associated



heterochromatin.,subunit:The nucleosome is a histone octamer containing two molecules each of H2A, H2B, H3 and H4 assembled in one H3-H4 heterotetramer

Background

Histones are basic nuclear proteins that are responsible for the nucleosome structure of the chromosomal fiber in eukaryotes. Nucleosomes consist of approximately 146 bp of DNA wrapped around a histone octamer composed of pairs of each of the four core histones (H2A, H2B, H3, and H4). The chromatin fiber is further compacted through the interaction of a linker histone, H1, with the DNA between the nucleosomes to form higher order chromatin structures. This gene encodes a replication-independent histone that is a member of the histone H2A family. It replaces conventional H2A histones in a subset of nucleosomes where it represses transcription and may participate in stable X chromosome inactivation. [provided by RefSeq, Oct 2015],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from 3T3 cells, primary antibody was diluted at 1:1000, 4° over night