




LRC4B rabbit pAb

Catalog No	YP-Ab-08644
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	LRRC4B LRIG4
Immunogen	Synthesized peptide derived from human LRC4B AA range: 22-72
Specificity	This antibody detects endogenous levels of LRC4B at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LRC4B; LRRC4B LRIG4
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane; Single-pass membrane protein. Cell junction, synapse, presynaptic cell membrane .
Tissue Specificity	
Function	similarity: Contains 1 Ig-like C2-type (immunoglobulin-like) domain., similarity: Contains 9 LRR (leucine-rich) repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MCF-7 cells,
primary antibody was diluted at 1:1000, 4° over night