




OAF rabbit pAb

Catalog No	YP-Ab-08650
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	OAF NS5ATP13TP2
Immunogen	Synthesized peptide derived from human OAF AA range: 221-271
Specificity	This antibody detects endogenous levels of OAF at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Catabolin; IL 1; IL 1 beta; IL1 beta; IL1- beta; OAF; OAF NS5ATP13TP2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	extracellular exosome,
Tissue Specificity	
Function	similarity: Belongs to the OAF family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from CACO2 cells,
primary antibody was diluted at 1:1000, 4° over night