



CE104 rabbit pAb

Catalog No	YP-Ab-08651
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	CEP104 KIAA0562
Immunogen	Synthesized peptide derived from human CE104 AA range: 415-465
Specificity	This antibody detects endogenous levels of CE104 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CE14 HUMAN;Centrosomal protein of 14 kDa;Cep14;GlyBP;KIAA562; LOC246295;RP1-286D6.4;CE104;CEP104 KIAA0562
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell projection, cilium . Cytoplasm, cytoskeleton, microtubule organizing center, centrosome, centriole . Cytoplasm, cytoskeleton, microtubule organizing center, centrosome . Cytoplasm, cytoskeleton, spindle pole. In interphase non-ciliated cells, localizes to the distal ends of both the mother and daughter centrioles. In ciliated cells, present at the distal end of the daughter centriole, but not on the mother centriole, and at the tip of primary cilium. Localization at the ciliary tip is also observed in motile cilia. Throughout S phase, associated with both mother and daughter centrioles in each centrosome. During metaphase and telophase, present at both spindle poles. .
Tissue Specificity	
Function	alternative products:Additional isoforms may exist,sequence caution:Contaminating sequence. Potential poly-A sequence.,similarity:Contains 2 HEAT repeats.,



Background

This gene encodes a centrosomal protein required for ciliogenesis and for ciliary tip structural integrity. The mammalian protein contains three amino-terminal hydrophobic domains, two glycosylation sites, four cysteine-rich motifs, and two regions with homology to the glutamate receptor ionotropic, NMDA 1 protein. During ciliogenesis, the encoded protein translocates from the distal tips of the centrioles to the tip of the elongating cilium. Knockdown of the protein in human retinal pigment cells results in severe defects in ciliogenesis with structural deformities at the ciliary tips. Allelic variants of this gene are associated with the autosomal-recessive disorder Joubert syndrome, which is characterized by a distinctive mid-hindbrain and cerebellar malformation, oculomotor apraxia, irregular breathing, developmental delay, and ataxia. [provided by RefSeq, Feb 2016],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from HeLa cells, primary antibody was diluted at 1:1000, 4° over night