




CDRT4 rabbit pAb

Catalog No	YP-Ab-08663
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	CDRT4
Immunogen	Synthesized peptide derived from human CDRT4 AA range: 6-56
Specificity	This antibody detects endogenous levels of CDRT4 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CDRT4
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Expressed in fetal skeletal muscle and kidney.
Function	tissue specificity:Expressed in fetal skeletal muscle and kidney.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MCF-7 cells,
primary antibody was diluted at 1:1000, 4° over night