




GPTC4 rabbit pAb

Catalog No	YP-Ab-08677
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	GPATCH4 GPATC4
Immunogen	Synthesized peptide derived from human GPTC4 AA range: 247-297
Specificity	This antibody detects endogenous levels of GPTC4 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GPTC4;GPATCH4 GPATC4
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	similarity:Contains 1 G-patch domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HCT116 cells,
primary antibody was diluted at 1:1000, 4° over night