




DAPL1 rabbit pAb

Catalog No	YP-Ab-08678
Isotype	IgG
Reactivity	Human; Mouse
Applications	IHC;IF
Gene Name	DAPL1 EEDA
Immunogen	Synthesized peptide derived from human DAPL1 AA range: 49-99
Specificity	This antibody detects endogenous levels of DAPL1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:50-200. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DAPL1;DAPL1 EEDA
Observed Band	
Calculated Molecular Weight	
Cell Pathway	integral component of membrane,
Tissue Specificity	Expressed in hair follicle (at protein level).
Function	function:May play a role in the early stages of epithelial differentiation or in apoptosis.,tissue specificity:Expressed in hair follicle (at protein level).,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Immunohistochemical analysis of paraffin-embedded human uterus. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).