




UFC1 rabbit pAb

Catalog No	YP-Ab-08688
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	UFC1 CGI-126 HSPC155
Immunogen	Synthesized peptide derived from human UFC1 AA range: 90-140
Specificity	This antibody detects endogenous levels of UFC1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HSPC155; UFC1; Ufm1 conjugating enzyme 1; UFC1 CGI-126 HSPC155
Observed Band	
Calculated Molecular Weight	
Cell Pathway	cytoplasm, extracellular exosome,
Tissue Specificity	
Function	function: E2-like enzyme which forms an intermediate with UFM1 via a thioester linkage., similarity: Belongs to the ubiquitin-conjugating enzyme family. UFC1 subfamily.,
Background	UFC1 is an E2-like conjugating enzyme for ubiquitin-fold modifier-1 (UFM1; MIM 610553) (Komatsu et al., 2004 [PubMed 15071506]). [supplied by OMIM, Mar 2008],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from 293T cells,
primary antibody was diluted at 1:1000, 4°over night