




YIPF4 rabbit pAb

Catalog No	YP-Ab-08704
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	YIPF4 Nbla11189
Immunogen	Synthesized peptide derived from human YIPF4 AA range: 59-109
Specificity	This antibody detects endogenous levels of YIPF4 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FinGER4; Nbla11189; Protein YIPF4; Yip1 domain family; member 4; YIP1 family member 4; YIPF4; YIPF4 Nbla11189
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Golgi apparatus, cis-Golgi network membrane ; Multi-pass membrane protein .
Tissue Specificity	Expressed in keratinocytes (at protein level).
Function	similarity: Belongs to the YIP1 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HpeG2 cells,
primary antibody was diluted at 1:1000, 4°over night