




LRIQ3 rabbit pAb

Catalog No	YP-Ab-08705
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	LRRIQ3 LRRC44
Immunogen	Synthesized peptide derived from human LRIQ3 AA range: 417-467
Specificity	This antibody detects endogenous levels of LRIQ3 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LRIQ3; LRRIQ3 LRRC44
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	similarity: Contains 1 IQ domain., similarity: Contains 2 LRR (leucine-rich) repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from K562 cells,
primary antibody was diluted at 1:1000, 4° over night