




2B1B rabbit pAb

Catalog No	YP-Ab-08706
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	HLA-DRB1
Immunogen	Synthesized peptide derived from human 2B1B AA range: 214-264
Specificity	This antibody detects endogenous levels of 2B1B at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Alcohol sulfotransferase;EC 2.8.2.2;HSST2;Hydroxysteroid sulfotransferase 2; ST2B1;ST2B1_HUMAN;Sulfotransferase 2B1;Sulfotransferase family cytosolic 2B member 1;SULT2B1;SULT2B1a;SULT2B1b;2B1B;HLA-DRB1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from K562 cells,
primary antibody was diluted at 1:1000, 4° over night