




SPA24 rabbit pAb

Catalog No	YP-Ab-08714
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	SPATA24
Immunogen	Synthesized peptide derived from human SPA24 AA range: 92-142
Specificity	This antibody detects endogenous levels of SPA24 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SPA24; SPATA24
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Nucleus, nucleolus . Nucleus, nucleoplasm . Associated with chromatin. .
Tissue Specificity	
Function	function: May play a role in cytoplasm movement and removal during spermiogenesis.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from A431 cells,
primary antibody was diluted at 1:1000, 4° over night