




PHLB3 rabbit pAb

Catalog No	YP-Ab-08715
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	PHLDB3
Immunogen	Synthesized peptide derived from human PHLB3 AA range: 568-618
Specificity	This antibody detects endogenous levels of PHLB3 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PHLB3;PHLDB3
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	sequence caution:Translated as Glu.,similarity:Contains 1 PH domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from SH-SY5Y cells,
primary antibody was diluted at 1:1000, 4° over night