




H90B3 rabbit pAb

Catalog No	YP-Ab-08719
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	HSP90AB3P HSP90BC
Immunogen	Synthesized peptide derived from human H90B3 AA range: 511-561
Specificity	This antibody detects endogenous levels of H90B3 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	H90B3;HSP90AB3P HSP90BC
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm .
Tissue Specificity	
Function	caution:Could be the product of a pseudogene.,function:Molecular chaperone., similarity:Belongs to the heat shock protein 90 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MDA-MB cells,
primary antibody was diluted at 1:1000, 4° over night