




SNG4 rabbit pAb

Catalog No	YP-Ab-08721
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	SYNGR4
Immunogen	Synthesized peptide derived from human SNG4 AA range: 37-87
Specificity	This antibody detects endogenous levels of SNG4 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SNG4;SYNGR4
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane; Multi-pass membrane protein.
Tissue Specificity	
Function	similarity:Belongs to the synaptogyrin family.,similarity:Contains 1 MARVEL domain.,
Background	This gene encodes an integral membrane protein. The gene belongs to the synaptogyrin gene family. Like other members of the family the protein contains four transmembrane regions. The exact function of this protein is unclear. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HuvEc cells,
primary antibody was diluted at 1:1000, 4°over night