




3BP5L rabbit pAb

Catalog No	YP-Ab-08723
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	SH3BP5L KIAA1720 UNQ2766/PRO7133
Immunogen	Synthesized peptide derived from human 3BP5L AA range: 87-137
Specificity	This antibody detects endogenous levels of 3BP5L at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	3BP5L;SH3BP5L KIAA1720 UNQ2766/PRO7133
Observed Band	
Calculated Molecular Weight	
Cell Pathway	intracellular,cytoplasm,
Tissue Specificity	
Function	similarity:Belongs to the SH3BP5 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from CACO2 cells,
primary antibody was diluted at 1:1000, 4° over night