




SPT6H rabbit pAb

Catalog No	YP-Ab-08726
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	SUPT6H KIAA0162 SPT6H
Immunogen	Synthesized peptide derived from human SPT6H AA range: 938-988
Specificity	This antibody detects endogenous levels of SPT6H at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	emb 5;Histone chaperone suppressor of Ty6;hSPT6;KIAA0162;SPT6;SPT6H;SUPT6H KIAA0162 SPT6H
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	Ubiquitously expressed.
Function	function:Acts to stimulate transcriptional elongation by RNA polymerase II., similarity:Belongs to the SPT6 family.,similarity:Contains 1 S1 motif domain., similarity:Contains 1 SH2 domain.,subunit:Interacts with RNA polymerase II and the DRB sensitivity-inducing factor complex (DSIF complex), which is composed of SUPT5H and SUPT4H1.,tissue specificity:Ubiquitously expressed.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from CACO2 cells, primary antibody was diluted at 1:1000, 4° over night